

**ARCONIC - Aluminium base, RM, 6000, wrought, Disc 64 mm x 25 mm**

Art. ID ALCOA-WR-6000  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Chromium (Cr)     | [7440-47-3] | 0,01i               | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,02i               | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,25                | %    |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 0,9                 | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,08i               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,05i               | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,6                 | %    |        |        |
|                  | Titanium (Ti)     | [7440-32-6] | 0,03i               | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 0,08i               | %    |        |        |