

ARCONIC - Aluminium base, RM, 6000, wrought, Disc 64 mm x 25 mm

Art. ID ALCOA-WZ-6000
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Boron (B) | [7440-42-8] | <0,005 | % | | |
| | Chromium (Cr) | [7440-47-3] | <0,0050 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,03i | % | | |
| | Iron (Fe) | [7439-89-6] | 0,15i | % | | |
| | Gallium (Ga) | [7440-55-3] | 0,02i | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,7 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,01i | % | | |
| | Nickel (Ni) | [7440-02-0] | <0,0050 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,4 | % | | |
| | Titanium (Ti) | [7440-32-6] | <0,0050 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,005i | % | | |
| | Zinc (Zn) | [7440-66-6] | <0,0050 | % | | |