

## **Custom Solution - ICP multi-element standard solution VIII 24 components in HNO3 (2% )**

Art. ID LM24-CUS-70614  
Unit 125 mL  
Deliverydetails No Dangerous Good UN1760 III E1

### Description

Conservation: 24 months (originally closed bottle), 12 month after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Barium (Ba)	[7440-39-3]	100	mg/L		
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Bismuth (Bi)	[7440-69-9]	100	mg/L		
	Calcium (Ca)	[7440-70-2]	100	mg/L		
	Cadmium (Cd)	[7440-43-9]	100	mg/L		
	Cobalt (Co)	[7440-48-4]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Copper (Cu)	[7440-50-8]	100	mg/L		
	Iron (Fe)	[7439-89-6]	100	mg/L		
	Gallium (Ga)	[7440-55-3]	100	mg/L		
	Potassium (K)	[7440-09-7]	100	mg/L		
	Lithium (Li)	[7439-93-2]	100	mg/L		
	Magnesium (Mg)	[7439-95-4]	100	mg/L		
	Manganese (Mn)	[7439-96-5]	100	mg/L		
	Sodium (Na)	[7440-23-5]	100	mg/L		
	Nickel (Ni)	[7440-02-0]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		
	Strontium (Sr)	[7440-24-6]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		
	Thallium (Tl)	[7440-28-0]	100	mg/L		
	Zinc (Zn)	[7440-66-6]	100	mg/L		