

### **Custom Solution - 8 components in HNO3 5%**

Art. ID LM24-CUS-70616  
Unit 100 ml  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### **Description**

Conservation: Conservation 12 months. (6 months after opening)

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Bismuth (Bi)      | [7440-69-9] | 20                  | mg/L |        |        |
|                  | Calcium (Ca)      | [7440-70-2] | 2000                | mg/L |        |        |
|                  | Potassium (K)     | [7440-09-7] | 200                 | mg/L |        |        |
|                  | Lithium (Li)      | [7439-93-2] | 200                 | mg/L |        |        |
|                  | Magnesium (Mg)    | [7439-95-4] | 200                 | mg/L |        |        |
|                  | Sodium (Na)       | [7440-23-5] | 200                 | mg/L |        |        |
|                  | Osmium (Os)       | [7440-04-2] | 20                  | mg/L |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 20                  | mg/L |        |        |