

**CuSn4Zn2PS, Cylinder, diameter 40 mm x 30 mm**

Art. ID BAM-M397a  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified value	Zinc (Zn)	[7440-66-6]	1,87 ± 0,06	%		
Certified value	Lead (Pb)	[7439-92-1]	0,227 ± 0,008	%		
Certified value	Tin (Sn)	[7440-31-5]	3,98 ± 0,1	%		
Certified value	Nickel (Ni)	[7440-02-0]	0,337 ± 0,007	%		
Certified value	Antimony (Sb)	[7440-36-0]	0,097 ± 0,004	%		
Certified value	Sulfur (S)	[7704-34-9]	0,45 ± 0,04	%		
Value for information	Arsenic (As)	[7440-38-2]	2,9 ± 0,2	mg/kg		
Value for information	Selenium (Se)	[7782-49-2]	< 1	mg/kg		
Certified value	Tellurium (Te)	[13494-80-9]	< 1	mg/kg		