

Super-High-Energy Charpy V-Notch Specimens (NIST-Verification, 8-mm Striker) - Please ask for availability!

| | |
|-----------------|-----------------------------------|
| Art. ID | NIST-2098 |
| Unit | set |
| Deliverydetails | No Dangerous Good /not restricted |

Description

This Standard Reference Material (SRM) is intended for the verification of Charpy machines equipped with an 8-mm striker, in accordance with the current ASTM Standard E23 and the current ISO Standard 148-2. SRM 2098 is made from 9310 alloy steel. The bars are finished to length, heat-treated, and machined in SRM specimen lots ranging in size from 1200 to 2000 specimens. Each specimen has a lot number and an identification number (three or four digits). A unit of SRM 2098 consists of a set of five specimens needed to perform a single verification.