

Bingham Concrete Mixture for Rheological Measurements

Art. ID	NIST-2497
Unit	kit for one batch
Deliverydetails	No Dangerous Good

Description

This Standard Reference Material (SRM) is intended for use in calibrating rheometers for measuring the rheological properties of concrete. A unit of SRM 2497 consists of two 18.9 L (5-gallons) containers labeled P and A. Container P includes Limestone (12 kg) and 4 L of corn syrup. Container A includes one plastic bag each of 1 mm beads (6 kg) and 10 mm beads (20 kg).