

## **Recombinant Human Insulin-like Growth Factor 1 (Frozen)**

Art. ID	NIST-2926
Unit	3 vials x 0.25 mL
Deliverydetails	No Dangerous Good

### Description

This Standard Reference Material (SRM) is primarily intended for use in calibrating mass spectroscopy-based procedures and devices for the determination of Insulin-like Growth Factor 1 (IGF-1) in human serum. It can also be used for value-assignment of control materials. A unit of SRM 2926 consists of three vials, each containing approximately 0.25 mL of a solution of recombinant human IGF-1