

Custom Solution - 10 compounds in in HNO3 5%

Art. ID LM24-CUS-70666
Unit 1 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (12 months after opening)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Titanium (Ti)	[7440-32-6]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Vanadium (V)	[7440-62-2]	100	mg/L		
	Phosphorus (P)	[7723-14-0]	100	mg/L		
	Boron (B)	[7440-42-8]	100	mg/L		
	Molybdenum (Mo)	[7439-98-7]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		
	Beryllium (Be)	[7440-41-7]	100	mg/L		