

Copper base, CuZn 39 Pb 3 , Setting-up-sample, disc Ø 40 mm x 25 mm

Art. ID SUS-2.0401
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Please note that the values varies from lot to lot.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zinc (Zn)	[7440-66-6]	39.32	%		
	Lead (Pb)	[7439-92-1]	3,21	%		
	Tin (Sn)	[7440-31-5]	0,19	%		
	Phosphorus (P)	[7723-14-0]	0,007	%		
	Manganese (Mn)	[7439-96-5]	0,009	%		
	Iron (Fe)	[7439-89-6]	0,26	%		
	Nickel (Ni)	[7440-02-0]	0,1	%		
	Arsenic (As)	[7440-38-2]	<0,01	%		
	Silver (Ag)	[7440-22-4]	0,008	%		
	Sulfur (S)	[7704-34-9]	0,002	%		