

Aluminium base, AlMg 3 Mn 0,3, Setting-up-sample, disc Ø 40 mm x 25 mm

Art. ID SUS-3.3525
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Please note that the values varies from lot to lot.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silicon (Si)	[7440-21-3]	0,08	%		
	Iron (Fe)	[7439-89-6]	0,17	%		
	Copper (Cu)	[7440-50-8]	0,004	%		
	Manganese (Mn)	[7439-96-5]	0,22	%		
	Magnesium (Mg)	[7439-95-4]	2,8	%		
	Chromium (Cr)	[7440-47-3]	0,001	%		
	Nickel (Ni)	[7440-02-0]	0,002	%		
	Zinc (Zn)	[7440-66-6]	0,007	%		
	Titanium (Ti)	[7440-32-6]	0,009	%		
	Lead (Pb)	[7439-92-1]	0,002	%		
	Tin (Sn)	[7440-31-5]	0,002	%		