

Custom Solution - Diuron 100 mg/mL in Methanol

Art. ID LM24-CUS-10132
Unit 2 x 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Haltbarkeit 18 Monat(e) (original verschlossene Flasche)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Diuron	[330-54-1]	100	µg/mL		