

### **Custom Solution - 2 components 1000 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-70694  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1000	mg/L		
	Barium (Ba)	[7440-39-3]	1000	mg/L		