

Custom Solution - 9 components 1000 mg/L each of in HNO3 5%

Art. ID LM24-CUS-70696
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

Description

Conservation 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	1000	mg/L		
	Cobalt (Co)	[7440-48-4]	1000	mg/L		
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Iron (Fe)	[7439-89-6]	1000	mg/L		
	Molybdenum (Mo)	[7439-98-7]	1000	mg/L		
	Nickel (Ni)	[7440-02-0]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		
	Thallium (Tl)	[7440-28-0]	1000	mg/L		
	Vanadium (V)	[7440-62-2]	1000	mg/L		