

### **Custom Solution - 5 components 100 mg/L each of in HNO3 2%**

Art. ID LM24-CUS-70697  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

#### Description

Conservation 24 months (originally closed bottle), 12 months after opening.

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Beryllium (Be)    | [7440-41-7] | 100                 | mg/L |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 100                 | mg/L |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 100                 | mg/L |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 100                 | mg/L |        |        |
|                  | Selenium (Se)     | [7782-49-2] | 100                 | mg/L |        |        |