

Custom Solution - 5 components 100 mg/L each of in HNO3 2%

Art. ID LM24-CUS-70697
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 24 months (originally closed bottle), 12 months after opening.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Chromium (Cr)	[7440-47-3]	100	mg/L		
	Lead (Pb)	[7439-92-1]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Selenium (Se)	[7782-49-2]	100	mg/L		