

## **High Chromium Cast Iron Block approx. 35 mm x 35 mm, thick.16 mm**

Art. ID DSZU-CH021  
Unit block  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,93 ± 0,06	%		
	Silicon (Si)	[7440-21-3]	0,52 ± 0,02	%		
	Mangan (Mn)		3,66 ± 0,09	%		
	Phosphorus (P)	[7723-14-0]	0,064 ± 0,003	%		
	Sulfur (S)	[7704-34-9]	0,009 ± 0,001	%		
	Chromium (Cr)	[7440-47-3]	9,07 ± 0,18	%		
	Nickel (Ni)	[7440-02-0]	5,86 ± 0,13	%		
	Titanium (Ti)	[7440-32-6]	0,093 ± 0,007	%		
	Copper (Cu)	[7440-50-8]	0,369 ± 0,017	%		
	Aluminium (Al)	[7429-90-5]	0,168 ± 0,029	%		
	Vanadium (V)	[7440-62-2]	0,61 ± 0,03	%		
	Molybdenum (Mo)	[7439-98-7]	4,42 ± 0,12	%		