

High Chromium Cast Iron Block approx. 35 mm x 35 mm, thick.18 mm

Art. ID DSZU-CH022
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,9 ± 0,06	%		
	Silicon (Si)	[7440-21-3]	0,43 ± 0,02	%		
	Mangan (Mn)		1,76 ± 0,05	%		
	Phosphorus (P)	[7723-14-0]	0,033 ± 0,002	%		
	Sulfur (S)	[7704-34-9]	0,018 ± 0,001	%		
	Chromium (Cr)	[7440-47-3]	14,85 ± 0,17	%		
	Nickel (Ni)	[7440-02-0]	2,19 ± 0,07	%		
	Titanium (Ti)	[7440-32-6]	0,078 ± 0,005	%		
	Copper (Cu)	[7440-50-8]	2,53 ± 0,05	%		
	Aluminium (Al)	[7429-90-5]	0,053 ± 0,015	%		
	Vanadium (V)	[7440-62-2]	0,45 ± 0,03	%		
	Molybdenum (Mo)	[7439-98-7]	2,65 ± 0,06	%		