

## **High Chromium Cast Iron Block approx. 35 mm x 35 mm, thick.17 mm**

Art. ID DSZU-CH024  
Unit block  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,01 ± 0,03	%		
	Silicon (Si)	[7440-21-3]	2,18 ± 0,06	%		
	Mangan (Mn)		1,22 ± 0,04	%		
	Phosphorus (P)	[7723-14-0]	0,102 ± 0,011	%		
	Sulfur (S)	[7704-34-9]	0,037 ± 0,004	%		
	Chromium (Cr)	[7440-47-3]	27,84 ± 0,36	%		
	Nickel (Ni)	[7440-02-0]	0,222 ± 0,009	%		
	Titanium (Ti)	[7440-32-6]	0,099 ± 0,009	%		
	Copper (Cu)	[7440-50-8]	0,88 ± 0,03	%		
	Aluminium (Al)	[7429-90-5]	0,096 ± 0,036	%		
	Vanadium (V)	[7440-62-2]	0,164 ± 0,009	%		
	Molybdenum (Mo)	[7439-98-7]	3,86 ± 0,22	%		