

## **High Chromium Cast Iron Block approx. 35 mm x 35 mm, thick.17 mm**

Art. ID DSZU-CH025  
Unit block  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	1,8 ± 0,05	%		
	Silicon (Si)	[7440-21-3]	2,7 ± 0,06	%		
	Mangan (Mn)		0,387 ± 0,021	%		
	Phosphorus (P)	[7723-14-0]	0,03 ± 0,004	%		
	Sulfur (S)	[7704-34-9]	0,026 ± 0,002	%		
	Chromium (Cr)	[7440-47-3]	35,14 ± 0,32	%		
	Nickel (Ni)	[7440-02-0]	1,77 ± 0,07	%		
	Titanium (Ti)	[7440-32-6]	0,117 ± 0,02	%		
	Copper (Cu)	[7440-50-8]	1,23 ± 0,03	%		
	Aluminium (Al)	[7429-90-5]	0,351 ± 0,028	%		
	Vanadium (V)	[7440-62-2]	0,044 ± 0,004	%		
	Molybdenum (Mo)	[7439-98-7]	0,302 ± 0,015	%		