

High Chromium Cast Iron Block approx. 35 mm x 35 mm, thick.18 mm

Art. ID DSZU-CH026
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	1,62 ± 0,04	%		
	Silicon (Si)	[7440-21-3]	1,14 ± 0,06	%		
	Mangan (Mn)		0,305 ± 0,018	%		
	Phosphorus (P)	[7723-14-0]	0,05 ± 0,003	%		
	Sulfur (S)	[7704-34-9]	0,032 ± 0,004	%		
	Chromium (Cr)	[7440-47-3]	35,87 ± 0,31	%		
	Nickel (Ni)	[7440-02-0]	3,63 ± 0,11	%		
	Titanium (Ti)	[7440-32-6]	0,013 ± 0,002	%		
	Copper (Cu)	[7440-50-8]	0,288 ± 0,019	%		
	Aluminium (Al)	[7429-90-5]	0,059 ± 0,018	%		
	Vanadium (V)	[7440-62-2]	0,067 ± 0,007	%		
	Molybdenum (Mo)	[7439-98-7]	0,96 ± 0,06	%		