

Iron Base, Cast Iron, RM, Mushroom 2 pieces 43 mm x 5 mm

Art. ID CTIF-GS5037
Unit 2 pieces
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	3,04	%		
	Chromium (Cr)	[7440-47-3]	0,014	%		
	Magnesium (Mg)	[7439-95-4]	0, ~0,184	%		
	Manganese (Mn)	[7439-96-5]	0,76	%		
	Nickel (Ni)	[7440-02-0]	0,64	%		
	Phosphorus (P)	[7723-14-0]	0,043	%		
	Sulfur (S)	[7704-34-9]	~0,0025	%		
	Silicon (Si)	[7440-21-3]	3,4	%		
	Titanium (Ti)	[7440-32-6]	0,029	%		