

Iron Base, Cast Iron, RM, Mushroom 2 pieces 43 mm x 5 mm

Art. ID CTIF-GS6135
Unit 2 pieces
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	~0,006	%		
	Carbon (C)	[7440-44-0]	3,6	%		
	Cobalt (Co)	[7440-48-4]	0,037	%		
	Chromium (Cr)	[7440-47-3]	0,04	%		
	Copper (Cu)	[7440-50-8]	0,0219	%		
	Magnesium (Mg)	[7439-95-4]	0, ~0,184	%		
	Manganese (Mn)	[7439-96-5]	0,38	%		
	Nickel (Ni)	[7440-02-0]	1,98	%		
	Phosphorus (P)	[7723-14-0]	0,013	%		
	Sulfur (S)	[7704-34-9]	~0,003	%		
	Silicon (Si)	[7440-21-3]	0,9	%		
	Titanium (Ti)	[7440-32-6]	0,007	%		
	Vanadium (V)	[7440-62-2]	0,0155	%		