

## **Iron Base, Cast Iron, RM, Mushroom 2 pieces 43 mm x 5 mm**

Art. ID CTIF-GS8018  
Unit 2 pieces  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,02	%		
	Carbon (C)	[7440-44-0]	3	%		
	Cerium (Ce)	[7440-45-1]	(<0,02)	%		
	Chromium (Cr)	[7440-47-3]	0,09	%		
	Copper (Cu)	[7440-50-8]	0,08	%		
	Magnesium (Mg)	[7439-95-4]	0, ~0,184	%		
	Manganese (Mn)	[7439-96-5]	0,7	%		
	Nickel (Ni)	[7440-02-0]	0,127	%		
	Phosphorus (P)	[7723-14-0]	0,07	%		
	Sulfur (S)	[7704-34-9]	~0,0015	%		
	Antimony (Sb)	[7440-36-0]	~0,01	%		
	Silicon (Si)	[7440-21-3]	3	%		
	Tin (Sn)	[7440-31-5]	0,07	%		
	Titanium (Ti)	[7440-32-6]	0,06	%		
	Vanadium (V)	[7440-62-2]	0,39	%		