

Iron Base, Cast Iron, RM , Mushroom 2 pieces 43 mm x 5 mm

Art. ID CTIF-FL2-1
Unit 2 pieces
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | ~0,006 | % | | |
| | Bismuth (Bi) | [7440-69-9] | ~0,0135 | % | | |
| | Carbon (C) | [7440-44-0] | 2,18 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,0263 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,44 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,0497 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,04 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | ~0,004 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,0238 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,049 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,082 | % | | |
| | Silicon (Si) | [7440-21-3] | 3,61 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,14 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,075 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,201 | % | | |