

## **Iron Base, Cast Iron, RM, Mushroom 2 pieces 43 mm x 5 mm**

Art. ID CTIF-FL4-1  
Unit 2 pieces  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	~0,05	%		
	Bismuth (Bi)	[7440-69-9]	~0,003	%		
	Carbon (C)	[7440-44-0]	2,6	%		
	Chromium (Cr)	[7440-47-3]	0,45	%		
	Copper (Cu)	[7440-50-8]	0,0168	%		
	Manganese (Mn)	[7439-96-5]	0,5	%		
	Molybdenum (Mo)	[7439-98-7]	0,09	%		
	Nitrogen (N)	[7727-37-9]	0,007	%		
	Nickel (Ni)	[7440-02-0]	0,061	%		
	Phosphorus (P)	[7723-14-0]	0,288	%		
	Sulfur (S)	[7704-34-9]	0,137	%		
	Silicon (Si)	[7440-21-3]	2,91	%		
	Tin (Sn)	[7440-31-5]	0,011	%		
	Titanium (Ti)	[7440-32-6]	0,0296	%		
	Vanadium (V)	[7440-62-2]	0,116	%		