

**Iron Base, Cast Iron, RM (FCr5), Mushroom 2 pieces 43 mm x 5 mm**

Art. ID CTIF-FCR5-1  
Unit 2 pieces  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Carbon (C)        | [7440-44-0] | 3,68                | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | 29,8                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 0,92                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,62                | %    |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7] | 3,15                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 2,49                | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,065               | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,027               | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,62                | %    |        |        |