

Copper base, Type: Bronze - Leaded, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-B11-1
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	84,75	%		
	Iron (Fe)	[7439-89-6]	0,17	%		
	Nickel (Ni)	[7440-02-0]	2	%		
	Phosphorus (P)	[7723-14-0]	0,057	%		
	Lead (Pb)	[7439-92-1]	1,93	%		
	Sulfur (S)	[7704-34-9]	0,09	%		
	Antimony (Sb)	[7440-36-0]	0,7	%		
	Silicon (Si)	[7440-21-3]	0,14	%		
	Tin (Sn)	[7440-31-5]	8,04	%		
	Zinc (Zn)	[7440-66-6]	2,1	%		