

Copper base, Type: Bronze - Tin, RM, 932 MOD, Mushroom 60 mm x 5 mm

Art. ID CTIF-B23-1
Unit each
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,02 | % | | |
| | Copper (Cu) | [7440-50-8] | 83,45 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,086 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,07 | % | | |
| | Lead (Pb) | [7439-92-1] | 7,2 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,019 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,384 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,025 | % | | |
| | Tin (Sn) | [7440-31-5] | 7,18 | % | | |
| | Zinc (Zn) | [7440-66-6] | 1,46 | % | | |