

**Copper base, Type: Bronze - Leaded, RM, Mushroom 60 mm x 5 mm**

Art. ID CTIF-B30-1  
Unit each  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,063	%		
	Copper (Cu)	[7440-50-8]	77,45	%		
	Iron (Fe)	[7439-89-6]	0,115	%		
	Manganese (Mn)	[7439-96-5]	0,15	%		
	Nickel (Ni)	[7440-02-0]	0,97	%		
	Phosphorus (P)	[7723-14-0]	0,063	%		
	Lead (Pb)	[7439-92-1]	10	%		
	Sulfur (S)	[7704-34-9]	0,048	%		
	Antimony (Sb)	[7440-36-0]	0,22	%		
	Silicon (Si)	[7440-21-3]	0,066	%		
	Tin (Sn)	[7440-31-5]	9,8	%		
	Zinc (Zn)	[7440-66-6]	1,05	%		