

**Copper base, Type: Brass, RM, Mushroom 60 mm x 5 mm**

Art. ID                      CTIF-L20-1  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,008	%		
	Arsenic (As)	[7440-38-2]	0,122	%		
	Copper (Cu)	[7440-50-8]	85,55	%		
	Iron (Fe)	[7439-89-6]	0,115	%		
	Manganese (Mn)	[7439-96-5]	0,043	%		
	Nickel (Ni)	[7440-02-0]	0,205	%		
	Lead (Pb)	[7439-92-1]	0,27	%		
	Silicon (Si)	[7440-21-3]	0,035	%		
	Tin (Sn)	[7440-31-5]	0,56	%		
	Zinc (Zn)	[7440-66-6]	13,1	%		