

Copper base, Type: Brass, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-UZ52-3
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]		%		
	Copper (Cu)	[7440-50-8]		%		
	Iron (Fe)	[7439-89-6]		%		
	Magnesium (Mg)	[7439-95-4]		%		
	Manganese (Mn)	[7439-96-5]		%		
	Nickel (Ni)	[7440-02-0]		%		
	Phosphorus (P)	[7723-14-0]		%		
	Lead (Pb)	[7439-92-1]		%		
	Antimony (Sb)	[7440-36-0]		%		
	Silicon (Si)	[7440-21-3]		%		
	Tin (Sn)	[7440-31-5]		%		
	Zinc (Zn)	[7440-66-6]		%		