

Iron Base, Cast Iron, RM, Mushroom 2 pieces 43 mm x 5 mm

Art. ID CTIF-FCR3-1
Unit 2 pieces
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,08 ± 0,023	%		
	Silicon (Si)	[7440-21-3]	0,25 ± 0,022	%		
	Manganese (Mn)	[7439-96-5]	1,02 ± 0,031	%		
	Sulfur (S)	[7704-34-9]	0,035 ± 0,002	%		
	Phosphorus (P)	[7723-14-0]	0,034 ± 0,003	%		
	Copper (Cu)	[7440-50-8]	0,0485 ± 0,002	%		
	Nickel (Ni)	[7440-02-0]	0,652 ± 0,007	%		
	Chromium (Cr)	[7440-47-3]	14,85 ± 0,14	%		
	Molybdenum (Mo)	[7439-98-7]	0,91 ± 0,02	%		