

Iron Base, Cast Iron, RM, Mushroom 2 discs Ø 43 mm x 5 mm

Art. ID CTIF-NR4-2S
Unit 2 discs
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Carbon (C) | [7440-44-0] | 2,47 | % | | |
| | Silicon (Si) | [7440-21-3] | 4,87 | % | | |
| | Manganese (Mn) | [7439-96-5] | 1,71 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,145 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,63 | % | | |
| | Nickel (Ni) | [7440-02-0] | 18,3 | % | | |
| | Chromium (Cr) | [7440-47-3] | 1,5 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,115 | % | | |