

## **Custom Solution - 23 components in HNO3 2% each at 10 mg/L**

Art. ID LM24-CUS-70723  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

### Description

Conservation 12 months. (6 months after opening).

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5]  | 10                  | mg/L |        |        |
|                  | Arsenic (As)      | [7440-38-2]  | 10                  | mg/L |        |        |
|                  | Boron (B)         | [7440-42-8]  | 10                  | mg/L |        |        |
|                  | Barium (Ba)       | [7440-39-3]  | 10                  | mg/L |        |        |
|                  | Beryllium (Be)    | [7440-41-7]  | 10                  | mg/L |        |        |
|                  | Cadmium (Cd)      | [7440-43-9]  | 10                  | mg/L |        |        |
|                  | Cobalt (Co)       | [7440-48-4]  | 10                  | mg/L |        |        |
|                  | Chromium (Cr)     | [7440-47-3]  | 10                  | mg/L |        |        |
|                  | Copper (Cu)       | [7440-50-8]  | 10                  | mg/L |        |        |
|                  | Iron (Fe)         | [7439-89-6]  | 10                  | mg/L |        |        |
|                  | Manganese (Mn)    | [7439-96-5]  | 10                  | mg/L |        |        |
|                  | Molybdenum (Mo)   | [7439-98-7]  | 10                  | mg/L |        |        |
|                  | Nickel (Ni)       | [7440-02-0]  | 10                  | mg/L |        |        |
|                  | Lead (Pb)         | [7439-92-1]  | 10                  | mg/L |        |        |
|                  | Antimony (Sb)     | [7440-36-0]  | 10                  | mg/L |        |        |
|                  | Selenium (Se)     | [7782-49-2]  | 10                  | mg/L |        |        |
|                  | Silicon (Si)      | [7440-21-3]  | 10                  | mg/L |        |        |
|                  | Strontium (Sr)    | [7440-24-6]  | 10                  | mg/L |        |        |
|                  | Tin (Sn)          | [7440-31-5]  | 10                  | mg/L |        |        |
|                  | Tellurium (Te)    | [13494-80-9] | 10                  | mg/L |        |        |
|                  | Thallium (Tl)     | [7440-28-0]  | 10                  | mg/L |        |        |
|                  | Vanadium (V)      | [7440-62-2]  | 10                  | mg/L |        |        |
|                  | Zinc (Zn)         | [7440-66-6]  | 10                  | mg/L |        |        |