

Custom Solution - 7 components 100 µg/mL each of in Ethylacetate

Art. ID LM24-CUS-70728
Unit 4 x 5 mL ampoule
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1173 II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-n-Octylphenol	[1806-26-4]	100	µg/mL		
	4-n-Octylphenol monoethoxylate	[51437-89-9]	100	µg/mL		
	4-n-Octylphenol-di-ethoxylate	[51437-90-2]	100	µg/mL		
	4-n-Nonylphenol	[104-40-5]	100	µg/mL		
	4-Nonylphenol-mono-ethoxylate	[104-35-8]	100	µg/mL		
	4-Nonylphenol-di-ethoxylate	[20427-84-3]	100	µg/mL		
	Bisphenol A D16	[96210-87-6]	100	µg/mL		