

### **Custom Solution - 3 components 1000 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-70745  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### Description

Conservation: 24 months (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	1000	mg/L		
	Lead (Pb)	[7439-92-1]	1000	mg/L		
	Tin (Sn)	[7440-31-5]	1000	mg/L		