

Custom Solution - 16 components, each 10 µg/mL in Methanol

Art. ID LM24-CUS-70762
Unit 3 x 1 mL (ampoule)
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Vinyl chloride	[75-01-4]	10	µg/mL		
	Chlorobenzene	[108-90-7]	10	µg/mL		
	Dichlorodifluoromethane	[75-71-8]	10	µg/mL		
	Fluorotrichloromethane	[75-69-4]	10	µg/mL		
	1,1-Dichloroethene	[75-35-4]	10	µg/mL		
	Methylene chloride (Dichloromethane)	[75-09-2]	10	µg/mL		
	trans-1,2-Dichloroethene	[156-60-5]	10	µg/mL		
	1,1-Dichloroethane	[75-34-3]	10	µg/mL		
	cis-1,2-Dichloroethene	[156-59-2]	10	µg/mL		
	1,2-Dichloropropane	[78-87-5]	10	µg/mL		
	Bromodichloromethane	[75-27-4]	10	µg/mL		
	Chlorodibromomethane	[124-48-1]	10	µg/mL		
	Bromoform (Tribromomethane)	[75-25-2]	10	µg/mL		
	1,3-Dichlorobenzene	[541-73-1]	10	µg/mL		
	1,2-Dichlorobenzene	[95-50-1]	10	µg/mL		
	1,2-Dibromo-3-chloropropane	[96-12-8]	10	µg/mL		