

### **Custom Solution - 3 components, each 1000 µg/mL in 2,2,4-Trimethylpentane (Isooctane)**

Art. ID LM24-CUS-70763  
Unit 3 x 1 mL (ampoule)  
Deliverydetails De Minimis UN1262 II E2

#### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Diisodecyl phthalate (mixture of isomers)	Diisodecyl phthalate (mixture of isomers)	[26761-40-0]	1000	µg/mL		
Diisononyl phthalate (mixture of branched chain isomers)	Di-iso-nonyl phthalate [28553-12-0]	[28553-12-0]	1000	µg/mL		
Phthalic acid, bis-isooctyl ester (technical)	Phthalic acid, bis-isooctyl ester (Bis-isooctyl phthalate)	[27554-26-3]	1000	µg/mL		