

Zinc Sulfide CRM

Art. ID AR-3085
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

Description

The intended use is for the Sulfur determination in materials using induction oxidation combustion with infrared detection or other valid test methods. The verification of this CRM was by induction oxidation combustion IR detection using tungsten metal and high purity iron chip for accelerators. This CRM was produced in accordance to ISO 17034 accreditation.

Lot #223F

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	32,9	%	Sulfur (S) by puriti y	
	Sulfur (S)	[7704-34-9]	32,9 ± 1,2	%	Sulfur (S) by analys is	