

## **Zinc Sulfide CRM**

Art. ID                      AR-3085  
Unit                         50 g  
Deliverydetails            No Dangerous Good /not restricted

### Description

The intended use is for the Sulfur determination in materials using induction oxidation combustion with infrared detection or other valid test methods. The verification of this CRM was by induction oxidation combustion IR detection using tungsten metal and high purity iron chip for accelerators. This CRM was produced in accordance to ISO 17034 accreditation.

Lot #223F

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur (S)	[7704-34-9]	32,9	%	Sulfur (S) by puriti y	
	Sulfur (S)	[7704-34-9]	32,9 ± 1,2	%	Sulfur (S) by analys is	