

**USP mAb 001, Monoclonal IgG1 - Biologic - (INTERNATIONAL COLD CHAIN SHIPMENT REQUIRED)**

Art. ID                      USP-1445539  
Unit                         each  
Deliverydetails            Dry ice shipment UN1845 E0

Description

This is a 200 µL solution of monoclonal antibody IgG1 in formulation buffer. Solution contains about 10.7 mg per mL of USP mAb 001 monoclonal IgG1 on the as is basis. Formulation buffer: 25 mM Citrate, 0.9% Sodium Chloride, 0.07% Polysorbate 80, pH 6.5

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	USP mAb 001, Monoclonal	[174722-31-7]				
	IgG1					