

**USP mAb 002, Monoclonal IgG1 - Biologic - (INTERNATIONAL COLD CHAIN SHIPMENT REQUIRED)**

Art. ID                      USP-1445547  
Unit                         each  
Deliverydetails            Dry ice shipment UN1845 E0

Description

This is a 200 µL solution of monoclonal antibody IgG1 in formulation buffer. Solution contains about 10.9 mg per mL of USP mAb 002 monoclonal IgG1 on the as is basis. Formulation buffer: 0.58% Monobasic sodium phosphate monohydrate, 0.12% Dibasic sodium phosphate anhydrous, 6% alpha,alpha-trehalose dihydrate , 0.04% Polysorbate 20, pH 6.2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	USP mAb 002, Monoclonal	[216974-75-3]				
	IgG1					