

**USP mAb 003, Monoclonal IgG1 - Biologic - (INTERNATIONAL COLD CHAIN SHIPMENT REQUIRED)**

Art. ID                    USP-1445595  
Unit                      each  
Deliverydetails        Dry ice shipment UN1845 E0

Description

This is a 200 µL solution of monoclonal antibody IgG1 in formulation buffer. Solution contains about 10.2 mg per mL of USP mAb 003 monoclonal IgG1 on the as is basis. Formulation buffer: 20 mM Histidine HCl, 120 mM Sucrose, 0.02% Polysorbate 20, pH 6.0 ± 0.3

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	USP mAb 003, Monoclonal	[912628-39-8]				
	IgG1					