

XRF Casting Dish (RFA Abgiessschale) Platinum (Pt) / Gold (Au) 95/5 % Ø 39/41mm x 3mm, Flange width 51,00mm, bottom thickness 0,80mm, Platinum 42,75 g / Gold 2,25 g

Art. ID	XRF-SM5
Unit	approx. 45 g
Deliverydetails	No Dangerous Good

Description

Used on the original Schoeps fusion machines as well as on machines of type HD-VAA, HD-Vulcan and Autofluxer
FAB - Casting Dish XRS Pt95/Au5 39/41x3mm SM5
RFA Abgiessschale aus Platin/ Gold 95/5 - XRS
passend für Autofluxer, Schoeps, Vulkan
Durchmesser 39/ 41 mm
Kantenlänge 51,00 mm
Innenhöhe 3,00 mm