

## **Custom Solution - 19 components in 5% HNO<sub>3</sub>**

Art. ID LM24-CUS-70803  
Unit 125 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

### Description

Conservation: 12 months (originally closed bottle).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Potassium (K)	[7440-09-7]	600	mg/L		
	Magnesium (Mg)	[7439-95-4]	600	mg/L		
	Boron (B)	[7440-42-8]	20	mg/L		
	Barium (Ba)	[7440-39-3]	20	mg/L		
	Copper (Cu)	[7440-50-8]	20	mg/L		
	Aluminium (Al)	[7429-90-5]	10	mg/L		
	Iron (Fe)	[7439-89-6]	10	mg/L		
	Zinc (Zn)	[7440-66-6]	10	mg/L		
	Manganese (Mn)	[7439-96-5]	4	mg/L		
	Arsenic (As)	[7440-38-2]	1	mg/L		
	Chromium (Cr)	[7440-47-3]	1	mg/L		
	Cadmium (Cd)	[7440-43-9]	1	mg/L		
	Nickel (Ni)	[7440-02-0]	1	mg/L		
	Lead (Pb)	[7439-92-1]	1	mg/L		
	Selenium (Se)	[7782-49-2]	1	mg/L		
	Antimony (Sb)	[7440-36-0]	1	mg/L		
	Thallium (Tl)	[7440-28-0]	1	mg/L		
	Uranium (U)	[7440-61-1]	1	mg/L		
	Vanadium (V)	[7440-62-2]	1	mg/L		