

405,9 µm Monodisperse Standards - NIST traceable

| | |
|-----------------|-----------------------------------|
| Art. ID | WHS-MS0406 |
| Unit | 5 bottle set, 1,5 g per bottle |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Ultra-narrow size distribution glass microspheres /// Only a few microns separating many grades /// Traceable to international standards - NIST & NPL /// The widest range of beads from a single source /// Over 50 sizes from 9 to 600µm /// Shape sorted for optimal sphericity /// Dry powders for excellent long term stability