

### **628 µm Monodisperse Standards - NIST traceable**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | WHS-MS0628                        |
| Unit            | 5 bottle set, 2,5 g per bottle    |
| Deliverydetails | No Dangerous Good /not restricted |

#### **Description**

Ultra-narrow size distribution glass microspheres /// Only a few microns separating many grades /// Traceable to international standards - NIST & NPL /// The widest range of beads from a single source /// Over 50 sizes from 9 to 600µm /// Shape sorted for optimal sphericity /// Dry powders for excellent long term stability