

Standard for 200 mm (8") diameter sieves - 800 µm, 850 µm, 900 µm, Mesh#20 (NIST Traceable)

| | |
|-----------------|-----------------------------------|
| Art. ID | WHS-SS412 |
| Unit | 5 bottle set, 2,5 g per bottle |
| Deliverydetails | No Dangerous Good /not restricted |

Description

Unique microsphere method of sieve calibration /// Mean aperture sizes traceable to NIST and NPL /// Measure over 1 million apertures per minute /// Typical calibration time of about 1 minute /// Method analyses over 80% of the sieve surface /// Accuracies and repeatabilities better than 1µm /// Results independent of sieve shaking method /// Single-shot bottles remove operator bias /// No need to send sieves away for calibration