

**Aniline - Mix 6, each je 100 µg/mL in Acetonitrile/Wasser (4:1)**

Art. ID LM24-CUS-10148  
Unit 1 x 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Conservation: 18 months (originally closed bottle)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	N-Methyl-2,4-dinitroaniline	[2044-88-4]	100	µg/mL		
	4-(2,4-Dinitrophenyl)morpholine	[39242-76-7]	100	µg/mL		
	N,N-dimethyl-2,4-dinitroaniline	[1670-17-3]	100	µg/mL		
	N-Ethyl-2,4-dinitroaniline	[3846-50-2]	100	µg/mL		
	2,4-Dinitrodiphenylamine	[961-68-2]	100	µg/mL		
	N-Cyclohexyl-2,4-dinitroaniline	[52790-66-6]	100	µg/mL		